

Compass Chair



It's the perfect inclusion chair. The rock-solid frame is built to withstand years of demanding use

The gracefully contoured seat and sturdy armrests provide reassuring boundaries, giving mildly involved students the support they need to feel secure and comfortable.

The five available sizes accommodate everyone from small children to adults, and the adjustable legs allow fine-tuning for optimal positioning.



The Compass Chair comes in five sizes.



DESIGNED FOR USE WITH THE
MOVE[®]
PROGRAM

Unique built-in support



The contoured sides of the seat, together with the inside surfaces of the armrests, provide subtle lateral support that can be particularly helpful for students with Autism Spectrum Disorder, giving security and sensory boundaries.



Optional Accessories



Seat belt


The optional seat belt clips into the seat when needed for safety or postural support.



Stability feet

The stability feet enlarge the chair's base of support. They are adjustable like the standard chair legs.

Compass Chair dimensions

User dimensions (cm)	R310 size 1	R320 size 2	R330 size 3	R340 size 4	R350 size 5
Height	91–119	109–137	124–152	140–173	157+
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.					
	Important: Make sure that seat width, depth, and height are adequate for each user, and that the user's weight does not exceed the maximum weight recommended.				
Item dimensions (cm)	R310	R320	R330	R340	R350
Seat height	23,25,28	27,29,32	32,34,37	37,39,42	42,44,47
Seat depth	24	28	32.5	37.5	42
Seat width	27	31	36	41.5	46.5
Backrest height	24	28	32	36.5	41
Armrest height above the seat	14	15	18	21	22.5
Armrest length	17	19	22	24.5	28
Max. working load (kg)	45	55	80	100	113

2019 order form

Use dimension chart to select appropriate size.

	R310 size 1	R320 size 2	R330 size 3	R340 size 4	R350 size 5
* Compass chair	R311 <input type="checkbox"/>	R321 <input type="checkbox"/>	R331 <input type="checkbox"/>	R341 <input type="checkbox"/>	R351 <input type="checkbox"/>
Seat belt	R312 <input type="checkbox"/>	R312 <input type="checkbox"/>	R332 <input type="checkbox"/>	R342 <input type="checkbox"/>	R342 <input type="checkbox"/>
Stability feet	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>
Total	R310	R320	R330	R340	R350

* Required

* Required

For all order enquiries please
 email admin@medix21.co.nz or
 Tel 04 233 1627

Exclusive to:
 Medix21
 29 Jarden Mile
 Ngauranga 6035
 Wellington

