

Hygienic, simple, affordable and adaptable to most toilets

Tool-free adjustments, assembly and disassembly make toileting and cleaning easier for caregivers.

Adaptable: Convert the HTS from stationary to mobile by removing the gray feet, and inserting casters.

Foam molded cushions: Optional soft foam pads for back and seat provide comfort for your clients and easy disinfecting for caregivers.

Hygienic design: The unique design of the seat pad opening – open to the rear – facilitates easy access for cleaning and hygiene.

Tilt-in-space: With its gas-assisted springs, the HTS tilts 5° forward for easier transfer, as well as 25° back for showering and hair washing.

Tilt-in-space: With its gas-assisted springs, the HTS tilts 5° forward for easier transfer,



DESIGNED FOR USE WITH THE  
**MOVE**  
PROGRAM

## Optional Accessories

Seat and Back Pads



Push handles



Headrest



Butterfly Harness



Anterior Support / Tray



Laterals



Laterals and Chest Strap



Hip Guides



Pan



Bowl Adapter



Deflectors



Abductor



Calf Rest and Ankle Straps



Footboard For use with Mounting Bar



Tub Base






# 2019 order form

Use dimension chart to select appropriate size.

		Z110 small	Z120 medium	Z130 large
<b>* Seat and back</b> (includes seat belt)		Z111 <input type="checkbox"/>	Z121 <input type="checkbox"/>	Z131 <input type="checkbox"/>
<b>* Base</b>	Mobile (tilt-in-space with footboard)	Z112 <input type="checkbox"/>	Z122 <input type="checkbox"/>	Z132 <input type="checkbox"/>
	Mobile (non-tilt with footboard)	Z113 <input type="checkbox"/>	Z123 <input type="checkbox"/>	Z133 <input type="checkbox"/>
	Stationary (non-tilt with footboard)	Z184 <input type="checkbox"/>	Z185 <input type="checkbox"/>	Z186 <input type="checkbox"/>
	Tub base (non-tilt)	Z138 <input type="checkbox"/>	Z138 <input type="checkbox"/>	N/A
	Mounting bar (round-front toilet)	Z115 <input type="checkbox"/>	Z115 <input type="checkbox"/>	Z135 <input type="checkbox"/>
	Mounting bar (elongated toilet)	Z125 <input type="checkbox"/>	Z125 <input type="checkbox"/>	Z135 <input type="checkbox"/>
<b>* Armrests</b> (pair)	Low	Z101 <input type="checkbox"/>	N/A	N/A
	Medium	Z102 <input type="checkbox"/>	Z102 <input type="checkbox"/>	Z102 <input type="checkbox"/>
	High	N/A	Z103 <input type="checkbox"/>	Z103 <input type="checkbox"/>
<b>Casters for mobile base</b>	Small casters (8 cm) special order only	Z195 <input type="checkbox"/>	Z195 <input type="checkbox"/>	Z195 <input type="checkbox"/>
	<b>New!</b> Large casters (10 cm)	Z181 <input type="checkbox"/>	Z181 <input type="checkbox"/>	Z181 <input type="checkbox"/>
<b>Seat and back pads</b> (recommended)	Rear-open seat pad	Z153 <input type="checkbox"/>	Z163 <input type="checkbox"/>	Z173 <input type="checkbox"/>
	Rear-closed seat pad	Z154 <input type="checkbox"/>	Z164 <input type="checkbox"/>	Z174 <input type="checkbox"/>
	No-hole seat pad	Z187 <input type="checkbox"/>	Z188 <input type="checkbox"/>	Z189 <input type="checkbox"/>
<b>Headrest</b>		Z104 <input type="checkbox"/>	Z104 <input type="checkbox"/>	Z104 <input type="checkbox"/>
<b>Push handles</b> (pair)		Z199 <input type="checkbox"/>	Z199 <input type="checkbox"/>	Z199 <input type="checkbox"/>
<b>Butterfly harness</b>		Z116 <input type="checkbox"/>	Z126 <input type="checkbox"/>	Z136 <input type="checkbox"/>
<b>Anterior support/tray</b>		Z145 <input type="checkbox"/>	Z155 <input type="checkbox"/>	Z165 <input type="checkbox"/>
<b>Laterals</b> (pair)	Small	Z105 <input type="checkbox"/>	Z105 <input type="checkbox"/>	Z105 <input type="checkbox"/>
	Large	N/A	Z106 <input type="checkbox"/>	Z106 <input type="checkbox"/>
	Small with strap	Z107 <input type="checkbox"/>	Z107 <input type="checkbox"/>	Z107 <input type="checkbox"/>
	Large with strap	N/A	Z108 <input type="checkbox"/>	Z108 <input type="checkbox"/>
<b>Hip guides</b> (pair)	Small	Z109 <input type="checkbox"/>	Z109 <input type="checkbox"/>	Z109 <input type="checkbox"/>
	Large	N/A	Z149 <input type="checkbox"/>	Z149 <input type="checkbox"/>
<b>Pan</b>		Z117 <input type="checkbox"/>	Z127 <input type="checkbox"/>	Z137 <input type="checkbox"/>
<b>Bowl adapter</b>		Z196 <input type="checkbox"/>	Z196 <input type="checkbox"/>	Z197 <input type="checkbox"/>
<b>Deflector</b>	Standard (requires seat pad)	Z118 <input type="checkbox"/>	Z118 <input type="checkbox"/>	Z128 <input type="checkbox"/>
	Splash guard	Z119 <input type="checkbox"/>	Z119 <input type="checkbox"/>	Z129 <input type="checkbox"/>
<b>Abductor</b>	Small	Z124 <input type="checkbox"/>	Z124 <input type="checkbox"/>	N/A
	Large	N/A	Z144 <input type="checkbox"/>	Z144 <input type="checkbox"/>
<b>Ankle straps</b> (pair, requires footboard)		Z156 <input type="checkbox"/>	Z156 <input type="checkbox"/>	Z157 <input type="checkbox"/>
<b>Calf rest</b> (requires footboard)		Z158 <input type="checkbox"/>	Z158 <input type="checkbox"/>	Z159 <input type="checkbox"/>
<b>Footboard</b> (for mounting bar attachment)		Z146 <input type="checkbox"/>	Z147 <input type="checkbox"/>	Z148 <input type="checkbox"/>
<b>Portability base with carry bag</b>		Z139 <input type="checkbox"/>	Z139 <input type="checkbox"/>	N/A
<b>Total</b>		<b>Z110</b>	<b>Z120</b>	<b>Z130</b>

\* Required

For all order enquiries please  
email  @medix21585 or  
Tel  W

p5.  
