

Hygienic, simple, affordable and adaptable to most toilets

Tool-free adjustments, assembly and disassembly make toileting and cleaning easier for caregivers.

Adaptable: Convert the HTS from stationary to mobile by removing the gray feet, and inserting casters.

Foam molded cushions: Optional soft foam pads for back and seat provide comfort for your clients and easy disinfecting for caregivers.

Hygienic design: The unique design of the seat pad opening – open to the rear – facilitates easy access for cleaning and hygiene.

Tilt-in-space: With its gas-assisted springs, the HTS tilts 5° forward for easier transfer, as well as 25° back for showering and hair washing.



DESIGNED FOR USE WITH THE
MOVE
PROGRAM

Optional Accessories

Seat and Back Pads



Push handles



Headrest



Butterfly Harness



Anterior Support / Tray



Laterals



Laterals and Chest Strap



Hip Guides



Pan



Bowl Adapter



Deflectors



Abductor



Calf Rest and Ankle Straps



Footboard
For use with
Mounting Bar



Tub Base



2019 order form

Use dimension chart to select appropriate size.

		Z110 small	Z120 medium	Z130 large
* Seat and back (includes seat belt)		Z111 <input type="checkbox"/>	Z121 <input type="checkbox"/>	Z131 <input type="checkbox"/>
* Base	Mobile (tilt-in-space with footboard)	Z112 <input type="checkbox"/>	Z122 <input type="checkbox"/>	Z132 <input type="checkbox"/>
	Mobile (non-tilt with footboard)	Z113 <input type="checkbox"/>	Z123 <input type="checkbox"/>	Z133 <input type="checkbox"/>
	Stationary (non-tilt with footboard)	Z184 <input type="checkbox"/>	Z185 <input type="checkbox"/>	Z186 <input type="checkbox"/>
	Tub base (non-tilt)	Z138 <input type="checkbox"/>	Z138 <input type="checkbox"/>	N/A
	Mounting bar (round-front toilet)	Z115 <input type="checkbox"/>	Z115 <input type="checkbox"/>	Z135 <input type="checkbox"/>
	Mounting bar (elongated toilet)	Z125 <input type="checkbox"/>	Z125 <input type="checkbox"/>	Z135 <input type="checkbox"/>
* Armrests (pair)	Low	Z101 <input type="checkbox"/>	N/A	N/A
	Medium	Z102 <input type="checkbox"/>	Z102 <input type="checkbox"/>	Z102 <input type="checkbox"/>
	High	N/A	Z103 <input type="checkbox"/>	Z103 <input type="checkbox"/>
Casters for mobile base	Small casters (8 cm) special order only	Z195 <input type="checkbox"/>	Z195 <input type="checkbox"/>	Z195 <input type="checkbox"/>
	New! Large casters (10 cm)	Z181 <input type="checkbox"/>	Z181 <input type="checkbox"/>	Z181 <input type="checkbox"/>
Seat and back pads (recommended)	Rear-open seat pad	Z153 <input type="checkbox"/>	Z163 <input type="checkbox"/>	Z173 <input type="checkbox"/>
	Rear-closed seat pad	Z154 <input type="checkbox"/>	Z164 <input type="checkbox"/>	Z174 <input type="checkbox"/>
	No-hole seat pad	Z187 <input type="checkbox"/>	Z188 <input type="checkbox"/>	Z189 <input type="checkbox"/>
Headrest		Z104 <input type="checkbox"/>	Z104 <input type="checkbox"/>	Z104 <input type="checkbox"/>
Push handles (pair)		Z199 <input type="checkbox"/>	Z199 <input type="checkbox"/>	Z199 <input type="checkbox"/>
Butterfly harness		Z116 <input type="checkbox"/>	Z126 <input type="checkbox"/>	Z136 <input type="checkbox"/>
Anterior support/tray		Z145 <input type="checkbox"/>	Z155 <input type="checkbox"/>	Z165 <input type="checkbox"/>
Laterals (pair)	Small	Z105 <input type="checkbox"/>	Z105 <input type="checkbox"/>	Z105 <input type="checkbox"/>
	Large	N/A	Z106 <input type="checkbox"/>	Z106 <input type="checkbox"/>
	Small with strap	Z107 <input type="checkbox"/>	Z107 <input type="checkbox"/>	Z107 <input type="checkbox"/>
	Large with strap	N/A	Z108 <input type="checkbox"/>	Z108 <input type="checkbox"/>
Hip guides (pair)	Small	Z109 <input type="checkbox"/>	Z109 <input type="checkbox"/>	Z109 <input type="checkbox"/>
	Large	N/A	Z149 <input type="checkbox"/>	Z149 <input type="checkbox"/>
Pan		Z117 <input type="checkbox"/>	Z127 <input type="checkbox"/>	Z137 <input type="checkbox"/>
Bowl adapter		Z196 <input type="checkbox"/>	Z196 <input type="checkbox"/>	Z197 <input type="checkbox"/>
Deflector	Standard (requires seat pad)	Z118 <input type="checkbox"/>	Z118 <input type="checkbox"/>	Z128 <input type="checkbox"/>
	Splash guard	Z119 <input type="checkbox"/>	Z119 <input type="checkbox"/>	Z129 <input type="checkbox"/>
Abductor	Small	Z124 <input type="checkbox"/>	Z124 <input type="checkbox"/>	N/A
	Large	N/A	Z144 <input type="checkbox"/>	Z144 <input type="checkbox"/>
Ankle straps (pair, requires footboard)		Z156 <input type="checkbox"/>	Z156 <input type="checkbox"/>	Z157 <input type="checkbox"/>
Calf rest (requires footboard)		Z158 <input type="checkbox"/>	Z158 <input type="checkbox"/>	Z159 <input type="checkbox"/>
Footboard (for mounting bar attachment)		Z146 <input type="checkbox"/>	Z147 <input type="checkbox"/>	Z148 <input type="checkbox"/>
Portability base with carry bag		Z139 <input type="checkbox"/>	Z139 <input type="checkbox"/>	N/A
Total		Z110	Z120	Z130

* Required

For all order enquiries please
email admin@medix21.co.nz or
Tel 04 233 1627

Exclusive to:
Medix21
29 Jarden Mile
Ngauranga 6035
Wellington

