

ECLIPSE



The *Eclipse* is a rigid chair designed for the bariatric user and **supports up to 1000 lbs** of occupant weight. Extended casters achieve **optimal weight balance** on the rear wheels for **excellent ease of propulsion**.



Hand Propulsion



Foot Propulsion



Self-Tilt



Desk/Table Accessible



Transit Tie Downs

Eclipse Rigid bariatric chair that turns on a dime

Available options See order form

- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ 30° manual dynamic recline
- ✓ Custom sizes (wider than 48")

visit pdgmobility.com/eclipse for more information



DYNAMIC TILT
No tilt



WEIGHT CAPACITY
Up to 1000 lbs



SEAT WIDTH
20" to 48"



SEAT DEPTH
16" to 24"



SEAT HEIGHT
14" to 20"



BACK SEAT ANGLE
85° to 100°



OVERALL WIDTH
+7.5"



CHAIR WEIGHT
44 lbs