

FUZE T50



The **Fuze T50** has a sleek aluminum frame with an infinitely adjustable tilt range of 50°. The chair is lightweight, has a **small footprint**, and is **easy to set up and adjust**. The T50 is PDG's best-selling tilt chair.



Hand
Propulsion



Foot
Propulsion



Self-Tilt



Desk/Table
Accessible



Transit Tie
Downs

Fuze T50 *The world's most versatile tilt chair*

visit pdgmobility.com/t50
for more information

The adaptable Fuze T50 offers a wide range of accessories and attachments. It is the world's most versatile tilt chair.



DYNAMIC TILT
Up to 50°



WEIGHT CAPACITY
Up to 350 lbs



SEAT WIDTH
12" to 24"



SEAT DEPTH
15" to 22"



SEAT HEIGHT
13" to 20"



BACK SEAT ANGLE
84° to 120°



OVERALL WIDTH
+4.0" to 8.5"



CHAIR WEIGHT
44 lbs



FUZE T50



Exclusive features

- ✓ A unique center of gravity pivot point
- ✓ Effortless and maintenance-free tilt performance
- ✓ Simple center of gravity adjustments

Available options *See order form*

- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ 30° manual dynamic recline
- ✓ Power tilt up to 350 lbs
- ✓ Ultra narrow overall width package
- ✓ Adjustable width frame



Ultra Narrow option 4" seat width plus





Available for
seat widths
18" to 24"

Some
limitations
apply



Power Tilt | 0 to full tilt in 20 seconds | up to 350 lbs



Hip Guides



Fold-down Backrest



Reclining Backrest | up to 30° of recline